



Anti-GST antibody DATA SHEET

Product Name: Mouse Anti-GST monoclonal antibody

Cat. No: BSABGST_010

Product type: primary antibodies

Description: Mouse monoclonal [6D] antibody to GST

Specificity: Recognizes GST recombinant proteins or GST proteins overexpressed in cells.

Immunogen: GST recombinant protein

Host: Mouse

Clonality: Monoclonal

Isotype: Mouse IgG1

Form: Liquid

Tested applications: WB, ICC/IF, ELISA, Lateral flow

Application: Western Blotting (WB, 1:5000-10000) · ELISA (1:5000-10000). Optimal working dilution for a specific application should be determined by the investigator.

Purity: Protein A purified

Components: Anti-GST purified antibody 0.1 mg

Storage buffer: Phosphate Buffered Saline, pH7.4

Storage instructions: 2-8°C