



Rabbit IgG polyclonal antibody DATA SHEET

Product Name: Goat Anti-Rabbit IgG polyclonal antibody

Cat. No: BS4002_100

Product type: primary antibodies

Description: Goat polyclonal antibody to Rabbit IgG

Specificity: BS4002_100 is specific for Rabbit IgG

Immunogen: Rabbit IgG, whole molecule

Host: Goat

Clonality: Polyclonal

Target: Rabbit IgG

Antigen Species: Rabbit

Species Reactivity: Rabbit

Tested applications: WB, ICC/IF, ELISA

Application: Western Blotting (WB, 1:5000-10000) · ELISA (1:5000-10000). Optimal working dilution for a specific application should be determined by the investigator.

Form: Purified polyclonal antibody from Goat serum with protein A column

Components: Anti-Rabbit IgG purified antibody 100 mg

Storage buffer: Phosphate Buffered Saline, pH7.4

Storage instructions: 2-8°C