



Rabbit IgG-HRP polyclonal antibody DATA SHEET

Product Name: Goat Anti-Rabbit IgG-HRP polyclonal antibody

Cat. No: BS4012_001

Product type: primary antibodies

Description: HRP-conjugated Goat anti-Rabbit IgG polyclonal antibody

Specificity: BS4012_001 is specific for Rabbit IgG-HRP

Immunogen: Rabbit IgG, whole molecule

Host: Goat

Clonality: Polyclonal

Target: Rabbit IgG

Antigen Species: Rabbit

Species Reactivity: Rabbit

Tested applications: WB, ICC/IF, ELISA

Application: Western Blotting (WB, 1:5000-10000) · ELISA (1:5000-10000). Optimal working dilution for a specific application should be determined by the investigator.

Conjugation: HRP

Form: Purified polyclonal antibody from Goat serum with protein A column

Components: Anti-Rabbit IgG-HRP purified antibody 1 mg

Storage buffer: 50 mM Tris-Glycine, pH7.4 , 0.15M NaCl, 40% Glycerol

Storage instructions: shipped at -20°C . Upon delivery aliquot and store at -20°C or -80°C . Avoid repeated freeze/thaw cycles.